



Dissolved Oxygen Meters

Jenway offer a bench dissolved oxygen meter with a wide range of accessories. In order to obtain additional descriptive information and specification details please select the appropriate product below.

Model	9500
Type	Bench
DO2 range % air saturation	0 to 199
DO2 range % oxygen saturation	0 to 25
DO2 range concentration (mg/l)	0 to 19.99
Temperature	Yes
Temperature units	°C/°F
Atmospheric pressure correction	Yes
Salinity correction	No
No. results stored	32
GLP compliant	No
Outputs	Analogue, RS232, IrDA