



## Vertical Autoclave (Steam Sterilizer)

Model:SH-AC-150M

Brand: SH Scientific

Origin: Korea

Feature:

- Automatic safety relief valve is installed for over pressure protection.
- Safe door opening mechanism & Plastic door cover for safety.
- Excellent temperature uniformity and accuracy.
- Stainless Steel Internal Chamber with baskets.
- Safety assurance with CE conformity.

## Specification

Model	SH-AC-150M
Temp Range	Max 130°C(266°F)
Max Pressure	0.12MPa, 15~16psi
Time Controller	Digital PID Controller
Time	00 to 99MIN (*Auto alert when the set time is done)
Pressure Gauge	Mechanical type 0 to 0.4 Mpa
Sensor	PT-100
Capacity	150L
Dimension Internal	510Øx740Hmm/ 120.1×29.1"
Dimension External (WxDxH)	670x690x1130mm / 26.4×27.2×44.5"
Net Weight	110Kg
Material Internal	Stainless Steel
Material External	Steel Plate with Powder Heating Coated
Heater	6KW
Electrical Requirements	27A
	220V, 50/60Hz, 1Ø