



Digital Water Baths,

Model: SWB6D

Brand: Stuart

Origin: England

Features:

- Choice of sizes: 6, 15 or 24 Liters
- Digital display and selection of temperature
- Low level water sensor
- Integral drain
- Supplied with robust polycarbonate lid

Three robust and reliable water baths, each with easy to clean stainless steel tank and sturdy metal outer casing. Digital water baths with easy to use temperature control and clear to read LED water temperature. Incorporated over temperature protection system that tracks the set temperature and controls the heater in the event of a fault. The heater is mounted underneath the tank to allow easy cleaning. Includes low level water sensor which cuts power to prevent the bath boiling dry. Incorporated drain for easy emptying of the bath. Minimum set temperature is 25°C.

Bath covers are included as standard with each bath to help reduce heat losses and evaporation. The covers are manufactured from robust polycarbonate allowing full visibility of the bath contents.

Specifications:

	SWB6D	SWB6D
Capacity		6L
Temperature range		+25°C to 99.9°C
Temperature stability		± 0.5°C
Internal dimensions (w x d x h),mm		300 x 150 x 200
Overall dimensions (w x d x h), mm		335 x 230 x 280
Net Weight, kg		6.5kg
Heater Power		350W
Electrical supply		230V, 50Hz

